




**Your trusted
global partner
in life-saving research**

For over 50 years, Allentown has provided critical solutions to the biomedical research community. Our high-quality equipment and advanced vivarium systems continue to define and exceed industry standards. Like everything we do, our unparalleled customer service is stamped with the Allentown hallmarks of integrity, dedication and care.

Sani-Rack™



AUTOCLAVE PRESENTATION RACK



OVERVIEW

Providing critical solutions with integrity and care.

Holds 28 cages:
4 wide x 7 high



Weight (unloaded):
148 lbs

Effective, Safe and Simple

The Sani-Rack is an effective, safe and simple method for autoclaving Sentry SPP or BCU-2 cage systems without the need for labor-intensive disassembly, stacking or costly cage bagging.



Stainless steel and high-temperature plastic construction

Silicone O-rings

Requires specific parameters during autoclaving

CERTIFIED HEPA cartridges for supply and exhaust ports (part # 230875-CS)

Autoclave cycle-count labels to track autoclave cycles

CERTIFIED HEPA cartridges can be used on as few, or as many cage slots as required, conserving HEPA life

Multiple cycles per HEPA*

Custom sizes available upon request