



NO FIT TESTING REQUIRED

MAXAIR Complete Systems are designed to help protect workers involved in critical research, biotechnology, and pharmaceutical manufacturing work when they may be exposed to potentially harmful airborne particulate contaminants.

---

**Our complete systems provide safer protection that is more comfortable to wear and easier to use and maintain.**

All MAXAIR SYSTEMS are NIOSH Approved



## SHROUD

Configured for situations where full 360° coverage of head and neck is desired and where reuse of filters protected by the Filter Cover Cap reduces on-going cost-per-use.

CA-DLC-CAPR-36<sup>1</sup> System with 2260-05ML/2260-05SM Single<sup>2</sup> Shroud



## CUFF

Comfortable and Convenient for moderate situations.

CA-DLC-CAPR-36<sup>1</sup> System with 2365-02ML/2365-02SM Cuff



## HOOD

Configured for situations calling for a single piece covering of the entire Helmet, Head, and Neck for maximum protection against self contamination and for fast donning and doffing.

CA-DLC-CAPR-36<sup>1</sup> System with 2271PB-07ML/2271PB-07SM Single<sup>2</sup> Hood

**\*MAXAIR Systems include numerous additional configurations to ensure there is a MAXAIR for your specific needs.**

<sup>1</sup>CA-DLC-DLC-37 12-15 Hr. Battery System optional in place of CA-DLC-CAPR-36 8-12 Hr System.

<sup>2</sup>Double Shroud 2261-01ML/2261-01SM and Double Hood 2272PB-07ML/2272PB-07SM available as options.